

QI GONG MACHINE SP

Dual outputs for stress reduction and pain management.

Alpha wave is known to promote body relaxation and natural healing.

Qi Gong Machine is based on the simulated human brain EEG and Chinese traditional medicine, Qi Gong (Chi Kung) (Book: Qi Gong For Health & Longevity, 1994, by Dr. Simon Wang and Dr. Julius Liu).

It can be applied to different parts of the body accordingly. The intensity can be selected at LOW or HIGH accordingly, such as comfort, body area and condition. LOW should be selected to the area of the head and the children, and HIGH for most applications. It is suggested to use to 10 or 30 minutes each session, 1-2 times per day for a period of 5-7 days.

This device is not to be substituted for the treatment of a physician. The safety of the device for use during pregnancy or birth has not been established. Pregnant women are urged to consult with their physician. Keep out of reach of children. This device is forbidden for use by the patient with a fever. Do not use while bathing or shower.

The device is not recommended for the patients with a heart disease. Do not place it over the pacemaker or over any implanted apparatus. Do not place it over the eyes.

The device can help on following conditions:

Acupressure

Apply the transducer to the acupuncture points

Intensity: Low or High

Time: 10 or 30 min/each time,
1-2 times/day.

Anxieties

Apply the transducer to the low abdomen area or to the center of soles (Acupuncture point K1)

Intensity: Low or High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Arthritis

Apply the transducer to the affected areas

Intensity: Low or High

Time: 10 or 30 min/each time,
1-2 times/day.

Acute Lumbar Sprain

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Back Pain

Apply the transducer to the back

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Carpal Tunnel Syndrome

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Chronic Fatigue Syndrome

Apply the transducer to the low abdomen area or to the center of soles (Acupuncture point K1)

Intensity: Low or High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Deep Tissue Massage

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day.

Energetic Bodywork

Apply the transducer to the low abdomen area or to the center of soles (Acupuncture point K1)

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Facial Paralysis

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Fibromyalgia

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Insomnia

Apply the transducer under the pillow or the mattress

Intensity: Low

Time: 10 or 30 min/each time.

Knee Injury

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Muscle Cramp

Apply the transducer to the affected areas

Intensity: Low or High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Muscle Sprain

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,
1-2 times/day for 5-8 days.

Osteoarthritis (hands and feet)

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Pain

Apply the transducer to the affected areas

Intensity: Low or High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Peripheral Neuropathy

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Post-stroke Rehabilitation

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Premenstrual Syndrome

Apply the transducer to the low abdomen area or the perineum or the center of soles (Acupuncture point K1)

Intensity: Low or High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Prostate Massage

Apply the transducer to the perineum

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Sciatica

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Soft Tissue Injury

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Sports Injury

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Sprain of the Ankle Joint

Apply the transducer to the affected areas

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Stiff Neck

Apply the transducer to the affected areas

Intensity: Low or High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Stress

Apply the transducer to the low abdomen area or the center of the soles (Acupuncture point K1)

Intensity: High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.

Tennis Elbow

Apply the transducer to the affected areas

Intensity: Low or High

Time: 10 or 30 min/each time,

1-2 times/day for 5-8 days.